



GLOVES

Revolution Maxi-Flex Gloves



Grey - 709 (displayed)

Now Available in ALL Revolution Colors!

683300-700	Black
683300-701	Red
683300-702	Fuchsia
683300-704	Pink
683300-705	Navy
683300-706	Teal
683300-707	Yellow
683300-708	Sky Blue
683300-709	Grey
683300-710	Purple
683300-711	White
683300-712	Ducks Green
683300-713	Tide Red
683300-714	Jewelry Box Blue
683300-715	Orange
683300-717	Light Blue Serenity
683300-718	Royal Blue
683300-719	Coral



Black - 700 (displayed)

SIZES:	15"
COLORS:	Revolution Colors
PROTECTION:	.50mm LE
ACCESSORIES:	683306 - Liners, 683307 - Covers

KEY FEATURES:

- Revolution fabric in solid colors
- .50mm LE
- Tough, and durable outer coverings
- Soft, flexible and comfortable
- Removable, washable and replaceable Cool Wear moisture wicking inner liners
- Replaceable outer covers

DESCRIPTION:

The perfect way to ensure your hands and lower forearms are fully protected by the most flexible, comfortable radiation protection gloves on the market. Comfortable fit with maximum dexterity. Soft removable, washable and replaceable inner liners reduce odors and allow for washing. Durable outer covers can be replaced. Available in 0.50 mm lead equivalency. Revolution Maxi-Flex 5 Finger Lead-Free Gloves come exclusively in our Revolution fabric, the softest and most comfortable material available for radiation protection.