



CURTAINS

Multiple Panel Table Mount Scatter Shield

077035



Universal adapter bracket fits almost every table with a side rail.

Model Number	Description
077035	Multiple Panel Scatter Shield

Upper Screen


- .50mm LE with stiff support
- 11" H x 26" W
- Hooks into mobile arm that can be easily removed

Lower Screen

- .50mm LE made up of 4 strips
- 30" L x 26" W

Additional Strip (Included)

- Applied with addition of intermediate joint
- 30" L x 10" W



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

