



BARRIERS

# Economy Mobile Curtain

## REV-MC



A unique new design for operating physicians that allows unparalleled protection for close-to-table studies. Sturdy, yet thoughtfully engineered with a small footprint so it can easily be positioned exactly where you want it for maximum scatter radiation protection.

- Raise and lower shield
- Smooth easy maneuvering
- .50mm LE double protective curtain
- 2 high-quality hospital grade locking casters

Model Number	Description
REV-MC	Economy Mobile Curtain
<b>Dimensions</b>	66" H x 28" W
<b>Curtain</b>	30" x 30" .50mm LE curtain with grey Revolution fabric
<b>Curtain Color</b>	Revolution Grey
<b>Frame Color</b>	White
<b>Frame</b>	Powder Coated Steel
<b>Wheels</b>	3 swivel casters, 2 with brakes
<b>Crate Size</b>	Box 1: 24" x 30" x 12" Box 2: 6" x 6" x 6"
<b>Crate Weight</b>	Box 1: 20 lbs / 9.07 kg Box 2: 27 lbs / 12.25 kg

