



BARRIERS

Large Window Mobile Barrier

683492



- Powder coated steel frames, precision TIG welded for years of trouble-free service completely enclosing the leaded glass panel for increased durability and safety
- Lead glass
- 1.5 mm lead equivalent opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.
- Shipped in white unless specified.

Model Number	Description
683492	Large Window Mobile Barrier - Lead
Dimensions	80¼" H x 40½" W X 25" D
Glass Size	72" H x 36" W
Lead Glass	2.0 mm Pb LE
Frame Color	White
Frame	Powder Coated Steel
Floor Curtain Color	Revolution Grey
Wheels	4 swivel casters, 2 with brakes
Crate Size	30" x 78" x 48"
Crate Weight	480 lbs / 217.72 kg