



BARRIERS

Regular Mobile Barrier

683484



- Powder coated steel frames, precision TIG welded for years of trouble free service completely enclosing the leaded glass panel for increased durability and safety
- Lead glass
- 1.5 mm lead equivalent opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.

Model Number	Description
683484	Regular Mobile Barrier - Lead
Dimensions	64½" H x 40½" W X 25" D
Glass Size	24" H x 36" W
Lead Glass	2.0 mm Pb LE
Frame Color	White
Frame	Powder Coated Steel
Curtain Color	Revolution Grey
Wheels	4 swivel casters, 2 with brakes
Crate Size	30" x 65" x 48"
Crate Weight	430 lbs / 195.04 kg



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

