



BARRIERS

Standard Mobile Barrier

683460



- Powder coated steel frames, precision TIG welded for years of trouble free service completely enclosing the leaded glass panel for increased durability and safety
- Lead glass
- 1.5 mm lead equivalent opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.

Model Number	Description
683460	Standard Mobile Barrier - Lead
Dimensions	72" H x 34½" W X 25" D
Glass Size	24" H x 30" W
Lead Glass	2.0 mm Pb LE
Frame Color	White
Frame	Powder Coated Steel
Floor Curtain Color	Revolution Grey
Wheels	4 swivel casters, 2 with brakes
Crate Size	30" x 72" x 42"
Crate Weight	335 lbs / 151.95 kg



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

