

Tilted Mobile Barrier

683461



- Powder coated steel frame
- Acrylic impact panels, enclosed in a precision TIG welded frame for years of trouble-free service, completely enclosing the leaded glass panel for increased durability and safety
- 2 high-quality hospital grade locking casters
- 1.50 mm lead equivalent opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.
- Ergonomical viewing area
- Shipping worldwide by ground, air or sea freight

Model Number	Description
683461	Tilted Mobile Barrier - Lead
683461-AC	Tilted Mobile Barrier - Acrylic
Dimensions	70" H x 34 1/4" W X 25" D
Glass Size	24" H x 30" W
Glass Options	Lead: 2.0 mm Pb LE Acrylic: .50 mm Pb LE
Frame Colors	White, Grey, or Black Custom colors available. Ships in standard white unless specified.
Frame Options	Powder Coated Steel Stainless Steel (Upgrade)
Curtain Color	Revolution Grey Custom colors available.
Wheels	4 swivel casters, 2 with brakes
Crate Size	73" x 30" x 43"
Crate Weight	355 lbs / 161.03 kg

