




BARRIERS

Revolution Mobile Hanging Acrylic Shield

REV-MHAS




Fully adjustable mobile shield that can be altered to user's preference. This product provides the benefits of an overhead suspension barrier with the mobility to fit into exam rooms of any size. The acrylic barrier can be raised, lowered, and rotated 360° in both vertical and horizontal planes to fit snug against the torso of the patient.

- Mobile unit with superior stability.
- Gives user full control of shield placement.
- Pivot ball allows for 360° motion.
- Durable 0.50mm acrylic shield.
- Small footprint, does not interfere with physician's position at the table.

Model Number	Description
REV-MHAS	Revolution Mobile Hanging Acrylic Shield

Dimensions	79" H x 37" W x 25" D
Window Size	24" H x 30" W
Protection	Acrylic window .50mm LE
Frame	Powder Coated Steel
Color	White
Wheels	4 swivel casters, 2 with brakes
Panel	.50mm LE
Crate Size	30" x 42" x 44"
Crate Weight	210 lbs / 95.25 kg



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

