

Dura-line™ XF+ flat panel detector

superior durability combined
with ultra-high resolution

Glass-free and high resolution Dura-line™ XF+ 10, XF+ 14, XF+ 17 wireless DR Detectors provide exquisite image quality, low weight and high durability.



High Image Quality at ultra-high resolution of 99 micron pixel

DQE 75% at 0 lp/mm



Superior durability

Glass-free substrate for
TFT drop height limit:
up to 1.5 meter



Convenient charging options

Cradle charger, USB-C charger,
tethered cable with magnetic
connector



IP67 Ingression protection

Up to 30 minutes safety operating
in water, and depth less than 1 meter



Greater portability

The ergonomic design
improves grips making it
easier to carry the detector
& preventing accidental drops



Stable wireless communication

IEEE802.11n/ac
2.4 GHz & 5 GHz



Long lasting batteries

Up to 15 hours operation
time (16 hours in stand by) -
when two batteries on board



Maximize your investment

Easy sharable between different
DR in-room and mobile modalities
and Retrofit, AED (Automatic
Exposure Detector) to support
vendorneutral DR Retrofit,
broad range of applications



Ergonomics

Low weight



Temperature stability

Operating temperature
0°C to 40°C



MUSICA®

Image processing for
excellent contrast detail and
examindependent, consistent
image quality



SmartXR™

AI assistant offers you
predictive workflow assistance,
for improved productivity and
consistent outcomes